Published online in: setox.eu/intertox & versita.com/science/medicine/it/ Copyright©2009 Slovak Toxicology Society SETOX



# Cholinesterase based amperometric biosensors for assay of anticholinergic compounds

#### **Miroslav POHANKA**

Center of Advanced Studies, Faculty of Military Health Sciences, University of Defence, Hradec Kralove, Czech Republic

ITX020209A05 • Received: 2 May 2009 • Revised: 10 May 2009 • Accepted: 15 May 2009

### ABSTRACT

Biosensors are analytical devices being approachable for multiple analytes assay. Here, biosensors with intercepted acetylcholinesterase (AChE) or butyrylcholinesterase (BChE) are presented as tool for assay of anticholinergic compounds such as pesticides, nerve agents and some natural toxins. Principle of assay is based on evaluation of cholinesterase activity and its pertinent decrease in presence of analyte. Nerve agents, pesticides, anticholinergic drugs useable for treatment of Alzheimer's disease as well as myasthenia gravis and aflatoxins are enlisted as compounds simply analyzable by cholinesterase biosensors.

**KEY WORDS:** biosensor; acetylcholinesterase; butyrylcholinesterase; amperometry; organophophates; organophosphonates; carbamates

### Introduction

Inhibition of two enzymes in body: acetylcholinesterase (AChE; EC 3.1.1.7) and butyrylcholinesterase (BChE; EC 3.1.1.8) is a common toxicology pathway of many compounds. Nerve agents, insecticides such as organophosphates as well as carbamates and some drugs may be named (Pohanka *et al.*, 2009a). The function of BChE is not crucial. On the other hand, AChE is an important enzyme in body. AChE is bound in outer membrane of neurons and terminates neurotransmission through neural cleft by hydrolyzing acetylcholine (Soreq and Seidman, 2001). Another action of AChE is the cholinergic anti-inflammatory pathway associated with vagus nerve (Pavlov *et al.*, 2003).

Intoxication by organophosphates, organophosphonates and carbamates is manifested by bronchospasm, bradycardia, miosis, lacrymation and diarrhoea. Furthermore, confusion, coma and/or respiratory failure caused by overstimulation of nicotinic and muscarinic receptors would be also observed when anticholinergic compound is administered (Eddieston *et al.*, 2008).

Correspondence address:

RNDr. Miroslav Pohanka, PhD.

Center of Advanced Studies, Faculty of Military Health Sciences, University of Defence, Trebesska 1575, 50002 Hradec Kralove, Czech Republic E-MAIL: rau@atlas.cz The inhibition of cholinesterases *in vitro* was found approachable for construction of analytical devices such as biosensors. The cholinesterase based biosensors are very effective tools for assaying of anticholinergic compounds. Even, estimation of *in vitro* impact of some drugs or presence of other anticholinergic compound is accompanied by inhibition of cholinesterases used in biosensor (Pohanka *et al.*, 2008a,b). The present review is aimed at summarization of the most important facts about cholinesterase based biosensors. Assay of typical anticholinergic compounds such as organophosphates, organophosphonates and carbamates by biosensors as well as biosensors principles and constructions are mentioned in the review.

### The most important cholinesterase inhibitors


There are available many different compounds with anticholinergic properties. Some of them, e.g. aflatoxins, are natural toxins with strong inhibition of AChE but not BChE (Egbunike and Ikegwuonu, 1984; Cometa et al., 2005; Pohanka et al., 2008e). Drugs (e.g. tacrine and galantamine being administered for treatment of Alzheimer's disease cognitive manifestation and myasthenia gravis are also modulators of AChE activity (Musial *et al.*, 2007). Some studies even found several plant secondary metabolites as strong inhibitors of cholinesterases (Lopez *et al.*, 2002).

Typical artificial inhibitors of cholinesterases are organophosphates and organophosphonates. These two

groups of anticholinergic compounds are irreversible inhibitors of both AChE and BChE. Toxic organophosphates are predominantly used as pesticides. Paraoxon, parathion, diisopropylfluorphosphate (DFP), malaoxon, malathion, chlorpyrifos and dichlorvos may be listed as typical examples (Marrs, 1993). Toxicity of some pesticides is attenuated by replacing of oxygen phosphate residuum by sulfur. It is e.g. derivate paraoxon and parathion or malaoxon and malathion. Toxicity is recovered due to mixed function oxidases in insect body. Warm blooded have not enzymatic apparatus for this pathway and thioderivates are slowly activated using cytochrome P-450 (Buratti et al., 2005). Some of organophosphonates were found useable as nerve agents in military points of view. Sarin, soman, tabun, cyclosarin and VX are well known representatives (Barthold et al., 2005). Examples of some above mentioned compounds are attached in Figure 1.

#### **Amperometric assay**

Amperometric biosensors are a good alternative to the potentiometric ones. Amperometric biosensors provide linear output signal dependence on concentration of analyte in comparison with the logarithmic relationship typical for potentiometric sensors (Ghindilis *et al.*, 1998). Amperometric biosensors are quite simple in comparison



with potentiometric devices. A simple noble metal wire can serve as the physicochemical transducer. Carbon paste electrodes are well known and also the screen-printed electrodes became one of the adopted possibilities. Carbon paste electrodes proved to have limited reproducibility due to coarsely defined carbon suspensions parameters. The best results were historically obtained with the screen-printed electrodes. The screen-printed electrodes can be simply produced on a mass scale. Miniaturization of electrodes seems to be ideal for flow-through based techniques (Pijanowska et al., 2003). Amperometric measurement can be realized in either two- or three-electrode configurations. The two electrode system is simpler. It consists of the reference and working electrodes. The two electrode configuration could be preferred for systems where currents or better current densities are low. If the value of current is high, the instability of reference electrode potential could occur and the three electrode configuration should be used. In this case, current flows between auxiliary and working electrodes.

The current going through the working electrode can be defined according to the Faraday law as follows:

$$i = \frac{Q}{t} = \frac{mFz}{Mwt}$$


The meaning of symbols is common; *i* means electric current, which is proportional to the electric charge Q transmitted during time interval *t*. The Faraday constant *F* is approximately equal to 96,485 C/mol. The transformed ions are described by molecular weight (*Mw*) and charge *z*; total mass of substances, which either precipitated or dissolved on the electrode, is expressed as *m*. An indifferent ion should be added to solution for polarization suppression.

## Amperometric evaluation of cholinesterase activity using biosensor

Biosensors are analytical devices consisting from biorecognition element and a proper sensor element (Brecht and Gauglitz, 1995). BChE and mainly AChE are promising recognition elements for biosensors construction (Pohanka *et al.*, 2008b). Evaluation of cholinesterase activity is the crucial factor in the construction of biosensors. Though the acetylcholine is commercially available, the reaction is not simply detectable. Typically pH electrodes would be applied as sensor detecting acidification of medium by releasing of acetic acid. Nowadays, glass electrodes are replaced by the semiconductors such as ion sensitive field effect transistor (ISFET) and light addressable potentiometric sensors (LAPS) being more approachable (Arkhypova *et al.*, 2001; Yoshinobu *et al.*, 2004). The replacement of potentiometric by amperometric principle has been found appropriate.

There is possibility to evaluate activity of cholinesterase by an amperometric principle in two ways. The first is based on performance of cholinesterase commonly with cholineoxidase (ChOx; EC 1.1.3.17) and oxygen or hydrogen peroxide amperometric transducer (Campanella *et al.*, 2007). ChOx oxidizes creating choline up betaine. A

# **54** Cholinesterase based amperometric biosensors for assay of anticholinergic compounds M. Pohanka



biosensor based on AChE and ChOx was found sensitive for assay of some pesticides such as pirimiphos-methyl in levels demanded by Europe Union (Del Carlo *et al.*, 2005).

The second way being frequently performed for amperometric evaluation of AChE activity is based on replacement of native substrate acetylcholine by an alternative acetylthiocholine. The mechanism is lucidly depicted in Scheme 1. Electrochemical oxidation of reaction product thiocholine is started by applied voltage (Pohanka *et al.*, 2009b). Biosensors with intercepted AChE working on above mentioned principle are promising devices for multiple assays. Pesticide paraoxon (Pohanka *et al.*, 2008c), dichlorvos (Valdes-Ramirez *et al.*, 2008), sulfotep (Kandimalla and Ju, 2006), natural toxic compound aflatoxin (Pohanka *et al.*, 2008e), nerve agents (Pohanka *et al.*, 2009b) and current as well as novel anticholinergic drugs (Pohanka *et al.*, 2008d) may be mentioned as examples of typical analytes.

### Conclusions

Amperometric biosensors with intercepted cholinesterases are promising tools for evaluation of many anticholinergic compounds. Assay of pesticides, nerve agents and aflatoxins are typical analytes. The second way of biosensor performance is evaluation of drugs administered to patients suffering from Alzheimer's disease and myasthenia gravis.

#### Acknowledgements

Ministry of Defence of Czech Republic is gratefully acknowledged for project No. OVUOFVZ200801.

#### REFERENCES

- Arkhypova VN, Dzyadevych SV, Soldatkin AP, Elskaya AV, Jaffrezic-Renault N, Jaffrezic H, Martelet C. (2001). Multibiosensor based on enzyme inhibition analysis for determination of different toxic substances. *Talanta* **55**: 919–927.
- Barthold CL and Schier JG. (2005). Organic phosphorus compounds nerve agents. *Crit Care Clin* **21**: 673–689.
- Brecht A and Gauglitz G. (1995). Optical probes and transducers. Biosens Bioelectron 10: 923–936.
- Buratti FM, D'Aniello A, Volpe MT, Menequz A and Testai E. (2005). Malathion bioactivation in the human liver: the contribution of different cytochrome p450 isoforms. *Drug Metab Dispos* **33**: 295–302.
- Campanella L, Lelo D, Martini E and Tomassetti M. (2007). Organophosphorus and carbamete pesticide analysis using an inhibition tyrosinase organic phase enzyme sensor; comparison by butyrylcholinesterase + choline oxidase opee and application to natural waters. *Anal Chim Acta* **587**: 22–32.

- Cometa MF, Lorenzini P, Fortuna S, Volpe MT, Meneguz A and Palmery M. (2005). In vitro inhibitory effect of aflatoxin B1 on acetylcholinesterase activity in mouse brain. *Toxicology* **206**: 125–135.
- Del Carlo M, Pepe A, Mascini M, De Gregorio M, Visconti A, Compagnone D. (2005). Determinating pirimiphos-methyl in durum wheat samples using an acetylcholinesterase inhibition assay. *Anal Bioanal Chem* **381**: 1367–1372.
- Eddieston M, Buckley NA, Eyer P and Dawson AH. (2008). Medical management of acute organophosphorus pesticide self-poisoning. *Lancet* **371**: 597–607.
- Egbunike GN and Ikegwuonu FI. (1984). Effect of aflatoxicosis on acetylcholinesterase activity in the brain and adenohypophysis of the male rat. *Neurosci Lett* **52**: 171–174.
- Ghindilis AL, Atanasov P, Wilkins M and Wilkins E. (1998). Immunosenors: electrochemical sensing and other engineering approaches. *Biosens Bioelectron* **13**: 113–131.
- Kandimalla VB and Ju H. (2006). Binding of acetylcholinesterase to multiwall carbon nanotube-cross-linked chitosan composite for flow-injection amperometric detection of an organophosphorous insecticide. *Chemistry* **12**: 1074–1080.
- Lopez S, Bastida J, Viladomat F and Codina C. (2002). Acetylcholinesterase inhibitory activity of some Amaryllidaceae alkaloids and Narcissus extracts. *Life Sci Narcissus* **71**: 2521–2529.
- Marrs TC. (1993). Organophosphate poisoning. Pharmacol Ther 58: 51-66.
- Musial A, Bajda M and Malawska B. (2007). Recent developments in cholinesterases inhibitors for Alzheimer's disease treatment. *Curr Med Chem* **14**: 2654– 2679.
- Pavlov VA, Wang H, Czura CJ, Friedman SG and Tracey KJ. (2003). The cholinergic anti-inflammatory pathway: a missing link in neuroimmunomodulation. *Mol Med* **9**: 125 – 134.
- Pijanowska DG, Sprenkels AJ, Olthuis W and Bergveld P. (2003). A flow-through amperometric sensor for micro-analytical systems. Sens Actuat B **91**: 98–102.
- Pohanka M, Jun D, Kalasz H and Kuca K. (2008b). Cholinesterase biosensor construction – a review. *Prot Pept Lett* **15**: 795–798.
- Pohanka M, Kuca K and Kassa J. (2008a). New performance of biosensor technology for Alzheimer's disease drugs: in vitro comparison of tacrine and 7-methoxytacrine. *Neuroendocrinol Lett* **29**: 755–758.
- Pohanka M, Musilek K and Kuca K. (2009a). Progress of biosensors based on cholinesterase inhibition. *Curr Med Chem* **16**: 1790–1798.
- Pohanka M, Dobes P, Drtinova L and Kuca K. (2009b). Nerve agents assay using cholinesterase based biosensor. *Electroanalysis* DOI: 10.1002/elan.200804528.
- Pohanka M, Jun D and Kuca K. (2008c). Amperometric biosensor for real time assay of organophosphates. *Sensors* 8: 5303–5312.
- Pohanka M, Kuca K and Kassa J. (2008d). New performance of biosensor technology for Alzheimer's disease drugs: in vitro comparison of tacrine and 7-methoxytacrine. *Neuroendocrinol Lett* **29**: 755–758.
- Pohanka M, Kuca K and Jun D. (2008e). Aflatoxin assay using an amperometric sensor strip and acetylcholinesterase as recognition element. *Sens Lett* **6**: 450–453.
- Soreq H and Seidman S. (2001). Acetylcholinesterase new role for an old actor. Natur Rev 2: 294 – 302.
- Valdes-Ramirez G, Fournier D, Ramirez-Silva MT and Marty JL. (2008). Sensitive amperometric biosensor for dichlorvos quantification: application to detection of residues on apple skin. *Talanta* 74: 741–746.
- Yoshinobu T, Schoning MJ, Finger F, Moritz W and Iwasaki H. (2004). Fabrication of thin film LAPS with amorphous silicon. *Sensors* **4**: 163–169.